# **Instructions Manual**

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# **CLINIC SW-S101**

**EN**

Please read before installation and before use

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produced by

DMEGA Co., Ltd.

Address : 314 Bongwha-Ro (38-7 Jugyul-Dong), Siheung-  
City, Gyunggi-Do, Korea (Postal Code : 15046)

## Warning & Caution

WARNING AND CAUTION INSTRUCTION, PLEASE READ BEFORE OPERATION. While reading your manual, please pay close attention to areas labeled: WARNING AND CAUTION.

The following symbols of warning will be found on the pump.

	Electric shock (Attached label on the product)
---	--

	Hot surface (Attached label on the product)
---	---

	Alternating current
---	---------------------

	CE mark
---	---------

	Warning, consult accompanying documents
---	---

	Ground protection
---	-------------------

	“ON” (power: connection from the mains)
---	---

	“OFF” (power: disconnection from the mains)
---	---

	Date of manufacture
---	---------------------

	Manufacturer
---	--------------

	Serial number
--	---------------



	Refer to instructions manual
--	------------------------------

	Correct Disposal of Electric and Electronic Equipment
---	---

## Index

Description	5	Maintenance	9
Unpacking	5	Disposal	10
Mounting	5	Troubleshooting Chart	11
Installation Location	5	Warranty	12
Installation	6		
Operation	8		

## Part list

- |                               |                                |
|-------------------------------|--------------------------------|
| 1. Vacuum air connection tube | 10. Filter housing             |
| 2. Air / Liquid Separator     | 11. Electric control unit      |
| 3. Motor Impellor Chamber     | 12. Capacitor                  |
| 4. Separation unit case       | 13. Cooling fan cover          |
| 5. Waste water outlet         | 14. Exhaust air outlet         |
| 6. Auto by-pass drain tubing  | 15. Power On/Off switch        |
| 7. Rubber foot                | 16. Overload protective switch |
| 8. Suction Inlet hole         |                                |
| 9. Filter cover               |                                |

## Description

This semi-wet suction system has been tested and found to comply with the standards for medical devices in Electromagnetic Compatibility Directive 93 / 42 / EEC

- EN 55011 : 2009 / A1 : 2010, EN 60601-1-2 : 2007 / AC :2010
- EN61000-3-3: 2008

These standards are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other nearby devices. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the space between the equipment.
- Connect the equipment to an outlet on a circuit different from that to which the other devices are connected.
- Consult the manufacturer or field service technician for assistance.

## Unpacking

Carefully remove the machine from the shipping case. Preserve all paperwork for future reference. If damage has occurred from shipment, claim must be filed with the carrier immediately. Preserve the shipping carton for inspection by the carrier. If you are required to communicate with your dealer be sure to include your order numbers for quick identification.

## Mounting

Rubber feet are attached to the product. Rubber feet are excellent for applications involving a semi-flexible surface such as bench top; they help to isolate noise and eliminate creeping. However, the horizontal position of the motor shaft should be preferred.

## Installation Location

The suction system should be positioned preferably in a clean, dry and well-ventilated area. Please be sure not to block the ventilation ports located on the motor housing. The unit should be placed where the surrounding temperature remains between 10°C

and 40°C (50°F and 104°F). Always check to insure the chosen location is protected from direct or indirect moisture contact.

## **IMPORTANT!**

This semi-wet suction system is Class 1 Equipment and shall to be installed by qualified personnel.

Never operate this product if it has a damaged power cable or plug. If it is not working properly, has been dropped, damaged or has fallen into water, please contact the appropriate service center for inspection and repair.

Keep the power cable away from heated surfaces. All electrical products generate heat. To avoid serious burns, NEVER touch suction motor during or immediately after operation. The temperature to suction motor can reach 120°C.

Never block any air exhaust valve or position the device on a soft surface where the openings may be blocked. The air exhaust valve is for ventilation of the motor inside the housing. Keep the air exhaust valve free of lint, dirt and other foreign objects.


Use only in well-ventilated areas. Be sure to properly identify intake and discharge ports before using machine.

Do not disassemble. Disassembly or attempted repairs if done incorrectly can create electrical shock hazard or product damage. Defer servicing to qualified service centers only.

Do not use this product in or near areas where it may fall or be pulled into water or other liquids. Do not reach for this product if it has fallen into liquid. Never operate this product outdoors in the rain or in a wet area.


## **INSTALLATION**

### **Suction Inlet Hose Connection**


- Remove the cover of suction inlet and push-in tubing connector into the inlet opening.
- Connect suction inlet hose (Ø30, spiral hose) at tubing connector and secure with electrical insulation tape.

### Waste Outlet Hose Connection

- Use provided spiral hose with 19 mm interior diameter.
- Connect hose to the outlet opening with clockwise direction.
- Please secure with hose clip or electrical insulation tape.
- Always position higher than drainage to prevent backflow problem.

### Exhaust Air Hose Connection

- Connect exhaust air outlet hose into the opening.
- Stretch hose clamp and push in air outlet hose and reposition hose clamp to tighten.
- Always place end of the air outlet hose out of building or by a window to prevent

noisy problems, hot temperature and bad odor.

- Always keep clean water and moisture from end of exhaust air hose. (Be cautious, there is a chance of incoming water like rain or vapor when the exhaust air hose is longer than 3 meters.)

## OPERATION

### Start Up

Check electric power connection and that all tubing connections are correct before starting up the suction system.

Check if there are any foreign material in the cooling fan part of the motor. Check any abnormal noise or vibration of the suction system. If this happens, refer to the troubleshooting guide in this manual. After addressing the causes of the noise and vibration, restart the system.

### Hose connection

- Suction inlet Ø30 mm (internal)
- Waste outlet Ø19 mm (internal)
- Exhaust air hose Ø38mm (internal)

Specifications		
Main Voltage (VAC)	230	230
Frequency (Hz)	50	60
Electric Power (kW)	0.4	0.5
Current Consumption (A)	2.7	3.2
Vacuum regulated (mbar)	150	180
Max. Liquid Inlet (LPM)	5	5
Suction motor rotation speed (rpm)	2900	3440
Number of treatment units (unit)	1	1
Noise level (dB)	55	60
Weight (kg)	17	17
Dimensions L x W x H (mm)	280 x 540 x 570	280 x 540 x 570

## MAINTENANCE

### Periodic maintenance


Regular maintenance should be performed by dentist or staff as follows:

- Every year:
  - » Check the motor status to see if it is abnormal. If abnormal bearing sound of the motor is detected, contact authorized dealer for inspection.
  - » Check the plastic hose or joint whether any leakage problem is occurring. If necessary, replace parts as needed.
- Every 5 ~ 7 years:
  - » Replace the motor bearing, and rubber parts like O-Ring and gasket.

### Cleaning Suction inlet filter

Cleaning procedures of filter component.

- Pull out filter component from machine in a clockwise direction of cover.
- Clean filter in running water to remove particles or debris. (It is recommended to

clean it at least twice per month to avoid reduction of suction power.)

- Reconnect filter. (Please be sure that filter fits into filter housing opening correctly.)
- When the clinic facility closes for an extended period of time like summer holiday, please leave filter removed. It helps keep the separator dry.

### Cleaning Suction Unit

After every treatment – To avoid bad odor and reduce the possibility of infection, the

operator should draw/suck one glass of water before the next treatment.

After finishing last treatment of the day – To avoid bad odor and reduce the possibility of infection, the operator should clean the suction unit with NON-FOAMING cleaning and disinfectant solution.

### **CAUTION!**

DO NOT use a foaming cleaning & disinfectant agent.  
It causes backflow of water to the suction motor and can damage it.

### **Disposal**

The built-in electronic circuit board and other electric components should be disposed as electric waste. Other metal components may be disposed of as metallic waste.

If the product is returned to the local dealer, all connections should be closed so that they are watertight. Uncontaminated plastic components of the product may be recycled.

## Troubleshooting Chart

This troubleshooting is to be carried out by qualified technicians only.

Problem	Reason	Remedy
Unit does not work	Mains supply is not connected Incorrect electrical connection or power source	Check mains power, fuses and reset circuit breaker if necessary
Unusual noise has become louder	Solid particle or debris is entering the suction motor	Dismantle unit and clean the impeller of the motor
Excessive vibration	Damaged impeller Motor and/or impeller are dirty	Dismantle unit and clean the impeller of the motor
Unit is very hot	Wrong wiring Low voltage Suction inlet is clogged Operating pressure or vacuum is too high	Check wiring Supply proper voltage Clean primary filter Install a relief valve and pressure or vacuum gauge
Suction is too weak	Suction inlet is clogged Leaking problem in the suction plumbing	Clean primary filter Check all suction tubes, hoses and connection, and replace if necessary
Suction power is weak and wastewater discharges during suction	Red rubber drain valve is in backside position or fitted incorrectly	Unscrew 2 bolts and position rubber drain valve correctly. (Convex surface is outside and flat surface is inside.)
Motor overload	Low voltage	Check power source Check wire size and wire connections

## WARRANTY

1. **WARRANTY DECLARATION:** the Manufacturer guarantees to the End User the replacement, or repair or intervention relative to all appliances, where the machine has non-conformities which hinder the correct use and operation, only where such non-conformities are due to the fault of the Manufacturer. The Manufacturer reserves the right to adopt the best solution in order to restore the appliance conformity within a reasonable amount of time. Warranty does not cover the parts of the machine subject to normal wear and/or deterioration, or defects resulting from inadequate transport.
2. **END USER:** End User is the natural or juridical person who purchases the machine to directly use it. An End User who buys the machine as a natural person for purposes that do not fall within the scope of commercial or professional activity is defined as “consumer”. The End User is considered as a “professional or company” when the purchase is made by a professional or company acting for purposes falling within the scope of his commercial or professional activity.
3. **TERMS:** The Manufacturer is liable to the End User when the lack of conformity occurs within the 36 months from the purchase date. Replacements or repairs under warranty will not extend the warranty period as originally envisaged, the expiry of which remains unchanged.
4. **VALIDITY:** the warranty is valid if the End User notifies the Manufacturer of any lack of conformity within 2 months from the date on which the non conformity was detected, but in any case within the final term stated here-above in clause 3. The above mentioned notification must necessarily be accompanied by a regular proof of purchase (till receipt or invoice).
5. **EXPIRY:** the warranty expires in case of inappropriate use of the machine, or in case of violation of the instructions of the “Instruction Manual” supplied by the Manufacturer, or in case the non conformity is due to inexperience of the user or to a fortuitous accident. Use of the machine after the discovery of defects will free the Manufacturer from any warranty obligation or liability.
6. **RESPONSIBILITY:** The Manufacturer is exempted from any responsibility deriving from any damage to the End User, deriving from the loss or reduction of production, resulting from any lack of conformity.
7. **COSTS AND FEES:** The End User is entitled to obtain the restoration of the conformity of the machine under warranty without paying labor or material costs. The shipping methods are to be agreed with the authorized dealer or with the Manufacturer. In the event that the non-conformity is not due to the Manufacturer’s fault or in case the warranty has expired, the Manufacturer can charge the End User both with the transport and repair costs that will be indicated by the authorized dealer or by the Manufacturer.

Prodotto da

DMEGA Co., Ltd.

Address : 314 Bongwha-Ro (38-7 Jugyul-Dong), Siheung-  
City, Gyunggi-Do, Korea (Postal Code : 15046)

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